# cold Documentation

Release 1.0.0

Doga Gursoy

Apr 17, 2023

# **USER GUIDE**

1	Credits	3
2	Index	5
3	Indices and tables	11

Warning: Still under development. It can work but expect some hiccups.

COLD is a software package for processing Laue diffraction images collected with coded-apertures.

### CHAPTER

# ONE

# CREDITS

#### CHAPTER

#### TWO

### INDEX

### 2.1 Getting started

#### 2.1.1 Before you begin

Coming soon.

#### 2.1.2 Installation

Coming soon.

#### 2.1.3 Setting up runtime environment

Coming soon.

#### 2.1.4 A "Hello world!" example

Coming soon.

### 2.2 General remarks

#### 2.2.1 Data processing workflow structure

Coming soon.

#### 2.2.2 Packing and unpacking

#### 2.2.3 Additional resources

Coming soon.

# 2.3 Runtime cofiguration options

#### 2.3.1 Data readout options

```
file:
    path: myfolder'
    range: [0, 100]
    threshold: 10
    frame: [0, 2048, 0, 2048]
    ext: 'h5'
    chunks: 8
    type: 'stacked'
    h5:
        key: '/entry1/data/data'
```

#### 2.3.2 Geometry options

#### Aperture

```
geo:
   mask:
       material: 'Au'
       path: 'codes/code-debruijn-2-8-000.npy'
       pad: 300
       bitsizes: [15, 7.5] # [mu]
        resolution: 0.5
        thickness: 7.5 # [mu]
        smoothness: 0 # [mu]
       widening: 2.5 # [mu]
        dist: 1.16 # [mu]
        tiltx: 5 #
        cenx: 290 #
        tilty: 18 #
        ceny: 100 #
        step: 1 # [mu]
        calibrate:
            dist: [1.15, 1.25, 0.01] # [mm]
            tiltx: [0, 8, 1] # -6
            cenx: [0, 0, 20] # -180
            tilty: [0, 23, 1] # 23
            ceny: [-100, 100, 20]
```

Detector

```
geo:
    detector:
        shape: [2048, 2048] # [pixels]
        size: [409.6, 409.6] # [mm]
        rot: [-1.20139958, -1.21416739, -1.21878591] # [radian]
        pos: [28.871, 2.786, 513.140] # [mm]
```

Source

geo:
 source:
 grid: [-1.365, -1.156, 0.001] # [mm]

#### 2.3.3 Algorithm options

Decoding

```
algo:
    pos:
        method: 'lsqr'
        init: 0
    sig:
        method: 'nnls'
        init:
            maxsize: 120 # [mu]
            avgsize: 20 # [mu]
            atol: 4
```

**Peak searching** 

Indexing

# 2.4 Aperture focusing

## 2.5 Peak searching

Coming soon.

# 2.6 Indexing

Coming soon.

## 2.7 General remarks

#### 2.7.1 Licensing

We use the BSD-3 licence that grants anyone almost unlimited freedom to do with the software as they please.

#### 2.7.2 Coding style

Readability counts! Try to follow PEP 8 style guide as much as possible, and more importantly use common sense and be consistent.

# 2.8 File formats

#### 2.8.1 APS 34-ID-E data format

The data is stored as HDF5 file format. Some of the important keys are given below:

# 2.9 Data structures

#### 2.9.1 Measurement data

```
import cold
path = 'myfolder/'
data = cold.read(path)
print (data.dtype, data.shape)
```

#### 2.9.2 Coded-apertures

import cold
path = 'myfolder/mymask.npy'
mask = cold.mask(path)
print (mask.dtype, mask.shape)

# 2.10 Adding new models

Coming soon.

# 2.11 Parallel execution

Coming soon.

### 2.12 Future plans

#### CHAPTER

# THREE

# **INDICES AND TABLES**

- genindex
- modindex
- search